

— Ξ 0 E 2 0 Р S Ζ  $\mathbf{z}$ -Г  $\mathbf{Z}$ 0 S 2 ۲Ì Ъ

2

Seat Colors Black Vinyl

Body Color <sup>Red</sup>

> The LiteRider PTC personal transport chairs offer users two great cash retail power chair choices! Both models feature a unique tubular frame design, a comfortable, stadium style seat and 18Ah batteries for an operating range up to 9.5 miles. What makes them different? While both models allow for the seat and battery pack to be easily removed, the GP162 takes disassembly a step further with a base that disassembles quickly and easily!



## liteRider **ptc**

Model Number	GP160	GP162
Weight Capacity	300 lbs.	300 lbs.
Drive Wheels	Rear	Rear
Maximum Speed	3.5 mph	3.5 mph
Operating Range <sup>1</sup>	9.5 miles	9.5 miles
Ground Clearance	2.5″	2 5"
Turning Radius	26″	26″
Type Batteries	2-18 AH	2-18 AH
Colors	Red	Red
Freewheel Mode	Yes	Yes
Electronic Speed Control	Yes	Yes
Electro-Mechanical Brakes	Yes	Yes
Charger	Off Board	Off Board
Controller	PG VR-2 50 Amp Remote	PG VR-2 50 Amp Remote
Dimensions:		
Length	33″	33"
Width	22.5″	22.5″
Height (ground to top of back of seat)	33" - 34"	33" - 34"
Height (ground to top of back of headrest)	N/A	N/A
Ground to Top of Seat	20.5" - 21.5"	20.5" - 21.5"
Number of Seat Height Positions	2	2
Size of Increments	1.0″	1.0″
Weight of Unit: (assembled)	114 lbs.	116 lbs.
Base	67 lbs.	N/A
Front Half	N/A.	34 lbs.
Rear Half	N/A	35 lbs.
Seat w/Arms	22 lbs.	22 lbs.
Batteries	25 lbs/pack	25 lbs/pack
Tires:	Flat Free	Flat Free
Drive Wheels	9″	9″
Casters	6″	6″
Anti-Tip Wheels	3″	3″
Standard Seat:	Stadium Style	Stadium Style
Slide Seat	N/A	N/A
Back Height - no headrest	16″	16″
Back Height with headrest	N/A	N/A
Width x Depth	17″ x 16″	17″ x 16″
Color	Black Vinyl	Black Vinyl
Footrest:	Yes	Yes
Height adjustable	Yes	Yes
Angle adjustable	Yes	Yes
Warranty:	<b>F</b> ' <b>)</b> (	F: )/
	Five Years	Five Years
Frame		
Drive Train Electronics <sup>2</sup>	13 Months 13 Months	13 Months 13 Months

Disassembles easily for transport!



1. Remove seat



3. GP162 model only. Grasp the frame and locking handle; lift the handle to unlock, and lift the frame up and off the drive train.

2. Pull out battery pack




Perfect for transport!

Model GP162 shown.

## Available Accessories:

- Cane Holder Tube
- Oxygen Tank Holder
- Cup Holder
- Quad Cane Holder
- Walker Holder
- Crutch Holder
- Safety FlagPower Chair Cover
- Touch-up Paint
- Deluxe Pack N' Go
- Seat Belt (Optional)

Some accessories require a mounting bracket or mounting clips.



<sup>2.</sup> Electronics warranty excludes batteries.

<sup>1.</sup> Battery range will vary due to rider weight, drive surface and terrain.



Golden Technologies • 401 Bridge Street, Old Forge PA 18518 • Tel: (800) 624-6374 • Fax: (800) 628-5165 • www.goldentech.com