

SHOWER BUDDY SB2

Tub Base (at bottom of legs) W x L: 465mm L x 335mm W (18.3" x 13.2")

PRODUCT SPECS WEIGHTS & DIMENSIONS

Weight Limit: 350 lb (160kg)

Rolling Base (W x D): 19.3" x 26.7" (490mm x 680mm)

Rolling Base including Footrest

19.3" x 30.7" (490mm x 780mm)

Outer Armrest Width: 22.5" (570mm)

Inner Armrest Width: 18.7" (475mm)

Inner Width Between Main Legs: 15.7" (400mm)

Backrest Width: 18" (460mm)

Back Rest Height from Cushion: 17" (440mm)

Toilet Height Clearance (Max): 19.4" (495mm)

Standard Bridge Length (Can Be Customized): 11.8" (300mm)

Seat Cushion WxD: 17.7" x 16.5" (450mm x 420mm)

Seat Cushion Height from Floor (Adjustable): 22.4" – 26.4" (570mm - 670mm)

Seat Cushion Opening (W x L): 6.8" x 13" (175mm x 335mm)

Tub Base (Top) WxL: 13.2" x 18.3" (335mm x 465mm)

Tub Base (side stabilizers): 16.7" – 28.3" (425mm – 720mm)

PRODUCT FEATURES

DETAILS

Armrests: Removable / Lockable / Flip-back

Height: Height Adjustable

Frame Metal: Corrosion resistant Aluminum frame (Grade 6061 T6)

Misc. Parts: Stainless steel parts and hardware (Grade 304)

Safety Belts: Lap & Chest Belts

Caster Wheels: 5" (125mm) , Lockable

Commode: Removable Commode System

SHIPPING INFORMATION SHIPPING WEIGHTS, DIMENSIONS

Weight: 67.00 LBS

Height: 31.00

Width: 22.00

Depth: 26.00

Sku: SB2

Warranty

WARRANTY

2 Year Limited Warranty: Showerbuddy warrants our Showerbuddy Solutions and all Showerbuddy products for two years from the date of purchase by the consumer. The warranty will only become active once the consumer submits the warranty registration card to the company with a copy of the purchase receipt. This warranty covers defects in materials and workmanship on all Showerbuddy products.

WARRANTY EXCLUSIONS

- Product or parts failure from normal wear and tear
- Damage caused by failure to adhere to the User's Manual
- Damage from improper installation, use or repair
- Brake line that operates the handbrake / tilt mechanism on tilting chairs (SB2T and SB3T)
- Improper use of the hand ratchet that is used in adjusting the neck support / adjustment
- Seat cushions and arm rests

LIMITATIONS

- The warranty is VOID if the original serial number tag is removed or altered.
- The warranty is VOID if the unit sustains damage from unnatural causes.
- The warranty is VOID if the unit is used to hold more than the weight limit of 300 pounds (140kg).
- The warranty is VOID if the unit has been subjected to neglect, improper maintenance, storage, cleaning or operation as defined in the User's Manual.
- The warranty is VOID if the unit is modified without Showerbuddy's written consent.
- This warranty applies to the original purchaser only, from the date of the purchase, provided a registration card has been submitted to the company within 60 days of purchase.

YOUR RESPONSIBILITIES

If your unit requires warranty service, please contact the dealer from which you purchased the Showerbuddy unit. If you do not receive satisfactory dealer service, contact Showerbuddy directly. Do not return the unit to Showerbuddy without first obtaining a Return Authorization Number from Customer Service. A defective unit must be shipped to Showerbuddy within 14 days of issuance of a Return Authorization Number by Customer Service.

SHOWERBUDDY'S RESPONSIBILITIES

Showerbuddy's sole responsibility is to repair or replace, at their discretion, the covered components of the unit.

This is the exclusive remedy for consequential or incidental damages. Upon pre-approval, Showerbuddy may provide exchange services for Warranty returns. Under this program, a refurbished Showerbuddy unit may be shipped after the return of the defective unit. The exchange unit may contain used parts that meet Showerbuddy quality standards. This program is provided at the sole discretion of Showerbuddy. The original warranty term will not be extended due to such a unit exchange and a return unit will become the property of Showerbuddy.

Returns

For hygienic reasons, this item is non-returnable.