


Brought to you locally by:

Call 1-888-637 8155 or visit us online to get started!

# Freecurve Dimensions

	SEAT	DIMENSIONS
A	Swivel radius	28 ¾" – 29 ½"
B	Wall to inner seat back	4 1⁄8" – 5 1⁄8"
С	Seat depth	14 ½" – 15"
D	Minimum folded width to footrest	16 ½" – 17 ¾"
E	Minimum open width to edge of footrest Offset footrest	25 1⁄8" – 26" 24" – 24 ¾"
F	Armrest width – external	22 1⁄8" – 23 1⁄2"
G	Armrest width – internal	17 1⁄8" - 18 1⁄4"
H	Seat back height	15 ½" – 16 ¾"
0	Footrest to seat height	18 5⁄8" – 19 ½"
J	Minimum footrest height	3 7%" / 4 7%"
	Minimum track intrusion into staircase	6 ¾" – 7 ½"
	Maximum weight limit	275 lbs

\* Dimensions vary based on seat style selected and may vary depending on hand of seat



Upholstery is made from a wipe clean and flame-retardant material.

All of our stairlifts are tested and approved by worldwide independent inspection services to give you the assurance that our products meet the highest standards in terms of quality and safety.





HC-BRO-0013



US Office 10888 Metro Court St. Louis, MO 63043 www.handicareusa.com Canadian Office 81 Romania Drive Concord, ON L4K 4Z9 www.handicare.ca Freecurve curved stairlift





# Freecurve

#### CURVED STAIRLIFT

A revolutionary design, the **Freecurve** utilizes a single tube that is custom built to fit your staircase. Coupled with a number of proprietary design options, the Freecurve maximizes your comfort while providing optimal aesthetics and convenience. Handicare has been a world leader in innovative lifting and mobility solutions for 130 years, providing safe solutions such as the Freecurve stairlift.

- LED illuminated footrest improves visibility in the dark, preventing you or other stair users from stumbling over the unfolded footrest
- Includes user-friendly retractable seatbelt
- Each unit comes with 2 handheld remotes to easily retrieve the stairlift
- Prevent unauthorized use with provided key
- Variety of seat options and track colors, including the ability to custom color match your track to your interior
- Lifetime of product, which is 10 years as stated on the manufacturers warranty documents.

### A tailor-made solution for your curved staircase

## Ergonomic and innovative design features

#### VERSATILE



A special seat for **small riders** ensure a safe & comfortable ride no matter their size.



Surrounded by safety sensors, Handicare stairlifts are designed to stop if meeting an obstruction - keeping you safe and still on vour seat.

#### **SUPPORTIVE**



With just a pull of the lever, the clever Active Seat option makes sitting down or standing up on the stairlift seat much easier by supporting your movements.

# Your first step to independence

#### SAFE

#### INTUITIVE



The easy-to-use joystick can be positioned on the left or right arm for **maximum control** and ensures a comfortable ride, every time.

#### ADAPTABLE



The Turn & Go feature automatically swivels your seat forward while riding. A convenient solution for staircases as narrow as 271/4" or users with limited knee flexion.



The Freecurve can be combined with a **powered hinge** or **drop nose** option to prevent the track from obstructing a doorway at the base of the stairs.



#### CONVENIENT