## Free-standing overhead systems

CANADA

PRISM MEDICAL

Waverley Glen free-standing overhead systems are everything free-standing systems were meant to be:

- they stand freely on the floor without support from floor-to-ceiling pressure posts
- they're collision resistant and don't require adjusting over time
- they're independent of ceiling type and wall type

The systems suit rental requirements as well as situations where ceiling or wall mounting is not possible or desirable. All free-standing systems are ideal for use with the Waverley Glen trans*portable* or Griffin lift. Units are assembled with standard wrenches and Allan keys.



## T-shape 3-post bedroom system

## Components

transfer rail with trolley
repositioning rail with trolley
cross braces
bed posts with leveling feet 1
room post with beveled foot

## Specifications

The system is 7 feet high without the leveling feet.

The repositioning track fits a bed length of up to 7 ft 3 in. The transfer track measures 4 ft 8 in.

Each post weighs 23.5 lbs. The tracks weigh 12 lbs and 8 lbs.